

Sanbonani Holdings (Pty) Ltd
RT2-1-2021 Supplementary Bid
Ending: 30 November 2024

michelle@sanbonani.co.za

emelia@sanbonani.co.za

Tel: 011 314 6617



RT2-05-36 - Domus 4
Alternating Pressure Redistribution System with Low Air Loss
R 15,442.63

Suitable for improving early-stage pressure injury prevention in long-term care patients.
Provide effective pressure relief care with alternating pressure redistribution and constant low-pressure modes.
Automatic adjustment for quick setup and pressure relief delivery.
Firmer air cells at the sacrum prevent patients from bottoming out when sitting.
Visual and audible alarm enhances patient safety and security.

5 Modes to Fulfil Individual Needs: Max Firm; Alternating; Continuous Low Pressure (CLP); Seat Inflation with Alternating; Seat Inflation with CLP

Mattress: 5" Overlay
Maximum Patient Weight: 200kg



Quick Connector
Transport Function



CPR Quick Release



Zipper-all-Around
Cover



RT2-05-37 - Domus Auto
Alternating Pressure Redistribution System With Automatic Adjustment
R 19,044.53

Domus Auto's automatic firmness calibration and manual refinement optimize early pressure injury prevention. Sacrum pressure relief allows individuals with different postures to sit longer and more comfortably.

5 Modes to Fulfil Individual Needs: Max Firm; Alternating; Continuous Low Pressure (CLP); Seat Inflation with Alternating; Seat Inflation with CLP

Mattress: 8" Replacement
Maximum Patient Weight: 250kg

RT2-05-38 - Pro-Care Auto 5"
Alternating Pressure Redistribution System
R 32,628.90

Suitable for patients in high to very high risk of pressure injuries.
 Alternating Low Pressure mode for an even higher comfort level to patients.
 Automatically provides an optimal pressure based on the patient's weight and position with a tuning function, Heel Relief and Transport Mode.
 Quiet Pump <20dBA
 Non-Slip PU Base with a Water Resistant and Flame retardant 4-way Stretch Cover.

Features: Auto Adjustment, Alternate Mode, Static Mode, Seat Inflate Mode, Max Firm Mode, Low Pressure Alarm, Power Failure Alarm, Alarm Mute, Panel Lock.

Mattress: 5" Overlay with a compartment for an additional foam mattress. (3")
 Maximum Patient Weight: 250kg



High to Very High Risk



Auto Sensing & Adjustment



Flame Retardant



Water Resistant



Fine-Tune Pressure




Visual & Audible Alarms



Vapour Permeable



Alternate Mode



RT2-05-39 - Pro-Care Auto 8"
Alternating Pressure Redistribution System
R 39,952.40

Suitable for patients in high to very high risk of pressure injuries.
 Alternating Low Pressure mode for an even higher comfort level to patients.
 Automatically provides an optimal pressure based on the patient's weight and position with a tuning function, Heel Relief and Transport Mode.
 Quiet Pump <20dBA
 Non-Slip PU Base with a Water Resistant and Flame retardant 4-way Stretch Cover.

Features: Auto Adjustment, Alternate Mode, Static Mode, Seat Inflate Mode, Max Firm Mode, Low Pressure Alarm, Power Failure Alarm, Alarm Mute, Panel Lock.

Mattress: 8" Replacement
 Maximum Patient Weight: 250kg

Sanbonani Holdings (Pty) Ltd
RT2-1-2021 Supplementary Bid
Ending: 30 November 2024



Sanbonani Holdings (Pty) Ltd

michelle@sanbonani.co.za
emelia@sanbonani.co.za

Tel: 011 314 6617