

Sanbonani Holdings (Pty) Ltd

Unit 4, 23 Venturi Crescent, Hennospark, Centurion, 0172

CSD: MAAA0005799

SAHPRA: 00000811MD

Company Registration No: 2004/000822/07

VAT Registration No: 4810216269

Practice No: 0195480



Enquiries:

Michelle Chamberlain

Tel: +27 11 314 6617

Fax: +27 11 314 6618

michelle@sanbonani.co.za

Specializing in Medical, Occupational Therapy and Physiotherapy Products



RT24-2021

SUPPLY, DELIVERY, INSTALLATION, COMMISSIONING AND
MAINTENANCE OF HOSPITAL FURNITURE AND WARD
REQUIREMENTS TO THE STATE FOR THE PERIOD 1 JUNE
2020 TO 31 MAY 2025

INDEX

Item Number	Description	Amount inclusive of 15% VAT	
RT24-08-001ME	Scale, Electronic, Infant – Cupid 1	R	1,304.05
RT24-09-001ME	Transfer Board, 80cm x 50cm – TMR6220 Medi-Roller Short with 50 disposable covers	R	4,207.63
RT24-09-002ME	Medicare Disposable Plastic Covers for TMR6220, pack of 50	R	641.66
RT24-09-003ME	Medicare Replacement Sliding Cover for TMR6220	R	725.82
RT24-09-004ME	Transfer Board, 180cm x 50cm – TMR6200	R	7,363.36
RT24-09-002ME	Medicare Disposable Plastic Covers for TMR6200, pack of 50	R	1,283.33
RT24-09-003ME	Medicare Replacement Sliding Cover for TMR6200	R	1,167.62

Please don't hesitate to contact us for more information or specifications.

Tel: 011 314 6617

michelle@sanbonani.co.za

emelia@sanbonani.co.za

lynette@sanbonani.co.za

sales@sanbonani.co.za

DIGITAL BABY SCALE – CUPID 1

RT24-08-001ME



CODE: SCAC1

Features:

- This balance is specially made for measuring the weight of a baby
- The pan may be removed to weigh babies under 20kg standing up
- This instrument is ideal for a family, hospitals or any infant care center
- Kg/lb switchable

Specifications:

MODEL:	CUPID 1
Capacity	0 – 20kg
	(0 - 10kg – 10g increments)
Overall Dimensions	320 (L) x 525 (W) x 160(H) mm
Display	1.0" LCD display
Platform Size	320 (L) x 520 (W) mm
Power Supply	9V Alkaline Battery AC Adapter
Unit Weight	2.2kg
Removable Tray	
Functions Keys	ON/OFF KG/LB, ZERO, HOLD, TARE
On/Off Button	Hold button down to turn unit off
Quality Standard	CE

Tel: 011 314 6617
Fax: 011 314 6618
www.sanbonani.co.za





Medi-roller Standard transfer board allows to move patients laterally comfortably and safely with less effort between beds, stretchers, operating tables, X-ray tables, etc.

Board made with an improved padded core more comfortable and light, covered with High Frequency welded polyamide fabric. Outer sliding sleeve made of very low friction silicone coated polyamide fabric. Also available disposable sleeve.

MEDI-ROLLER STANDARD

FEATURES



X-ray translucent and MRI safe.




Ultra sliding sleeve cover



High frequency welded board



Lightweight PE foam core



2 models: Large (foldable) Short

SPECIFICATIONS

PRODUCT FAMILY: Medi-roller

PRODUCT NAME: Medi-roller Standard

SIZE CHART

Ref.	Model	Measures (cm)	Weight (Kg)
TMR6200	Foldable	180 x 50	2,00
TMR6220	Short	90 x 50	1,00

COLOUR

- Blue board
- Orange cover

COMPOSITION

Board: Ext: 70% PVC - 30% Polyester
 Int: Polyethylene
Cover: Silicone coated Polyamide
LATEX FREE

LAUNDRY INSTRUCTIONS

Board: Clean with damp cloth and mild soap or surface disinfectant detergent.

Cover: - Machine wash max. 60°C - Do not tumble dry
 - Do not bleach - Do not iron
 - Do not dry clean




DISINFECTION RECOMMENDATIONS

- For infection prevention and control use the disposable protective cover (ref. TMR6270).
- Use a surface disinfectant or disinfectant wipes if required. For cover disinfection machine wash at 60°C using a disinfectant detergent. Do not use softener, chlorine nor whitening concentrates.

WARNING

- Maximum weight: 170 Kg.
- Maximum distance between surfaces: 20cm

PACKAGING



- 1 - Medi-roller Standard board
- 1 - Medi-roller Standard cover



Medi-roller Standard transfer board allows to move patients laterally comfortably and safely with less effort between beds, stretchers, operating tables, X-ray tables, etc.

Board made with an improved padded core more comfortable and light, covered with High Frequency welded polyamide fabric. Outer sliding sleeve made of very low friction silicone coated polyamide fabric. Also available disposable sleeve.

MEDI-ROLLER STANDARD

FEATURES



X-ray translucent and MRI safe.




Ultra sliding sleeve cover



High frequency welded board



Lightweight PE foam core



2 models: Large (foldable) Short

SPECIFICATIONS

PRODUCT FAMILY: Medi-roller

PRODUCT NAME: Medi-roller Standard

SIZE CHART

Ref.	Model	Measures (cm)	Weight (Kg)
TMR6200	Foldable	180 x 50	2,00
TMR6220	Short	90 x 50	1,00

COLOUR

- Blue board
- Orange cover

COMPOSITION

Board: Ext: 70% PVC - 30% Polyester
 Int: Polyethylene
Cover: Silicone coated Polyamide
LATEX FREE

LAUNDRY INSTRUCTIONS

Board: Clean with damp cloth and mild soap or surface disinfectant detergent.

Cover: - Machine wash max. 60°C - Do not tumble dry
 - Do not bleach - Do not iron
 - Do not dry clean




DISINFECTION RECOMMENDATIONS

- For infection prevention and control use the disposable protective cover (ref. TMR6270).
- Use a surface disinfectant or disinfectant wipes if required. For cover disinfection machine wash at 60°C using a disinfectant detergent. Do not use softener, chlorine nor whitening concentrates.

WARNING

- Maximum weight: 170 Kg.
- Maximum distance between surfaces: 20cm

PACKAGING



- 1 - Medi-roller Standard board
- 1 - Medi-roller Standard cover